




XENOSEP

FunnelFit™ Filters for Extract Drying Applications

Pre-Formed Conical Filter Saves Time and Solvent—Reduces Cost and Increases Lab Productivity

For about the same price as your current filter, new FunnelFit Conical Filters eliminate the mundane and frustrating preparation steps required when working with cumbersome, ill-fitting, pre-pleated or quadrant folded filters. Our pre-formed conical filter fits the funnel perfectly every time which improves quality by preventing sample bypass and tip blow-out.



DROP & GO

Using FunnelFit™ Conical Filters is 10 times faster than quadrant folding circular filters which improves laboratory productivity



FunnelFit™ Conical Filters fit the funnel perfectly—ensures high quality results



Eliminate folding, repositioning, and pre-wetting hassles to save time and solvent costs

Ordering Information

Catalog No.	Description/Grade	Qty/Pack	Compare to Whatman #
2440-125	125 mm DA-Medium (8-12 µm)	100	1440-125
2440-150	150 mm DA-Medium (8-12 µm)	100	1440-150
2441-125	125 mm DA-Coarse (20-25 µm)	100	1441-125
2441-150	150 mm DA-Coarse (20-25 µm)	100	1441-150
2442-125	125 mm DA-Fine (2-3 µm)	100	1442-125
2442-150	150 mm DA-Fine (2-3 µm)	100	1442-150

Other diameters and grades available upon request.

To obtain **FREE** samples and a competitive price quote contact...



Xenosep Technologies LLC
508 Main Street, Boonton, NJ 07005
866.XENOSEP (866.936.6737) • 973.239.0247 • FAX: 973.239.4133
E-mail: info@xenosep.com • www.xenosep.com

©Copyright, Xenosep Technologies, 2008 Printed in USA